

*Sample preparation for Permethylated GSL profiling
by MALDI/TOF-MS and NSI-MSn*

Date written: 10/8/18

Authors: Kazuhiro Aoki* and Mayumi Ishihara

Complex Carbohydrate Research Center, University of Georgia, Athens, GA 30602

*To whom correspondence should be addressed; Kazuhiro Aoki (kaoki@ccrc.uga.edu)

Reagents:

- Methanol (MeOH)
- α -dihydroxybenzoic acid (DHBA)
- HPLC Water
- Sodium Hydroxide Solution, 50% w/w
- Iso-propanol (2-PrOH)
- n-propanol (1-PrOH)

Equipment/Apparatuses:

- Glass culture tubes (round bottom)
- 500 μ L Microcentrifuge tube or glass inserts for GC autosampler vial
- GC glass autosampler vial (conical bottom)
- 1.5 mL Microcentrifuge tube
- Glass pipette
- Plastic pipette
- Glass syringe
- Centrifuge
- Vortex



Procedures:

MALDI/TOF-MS analysis

MALDI/TOF-MS analysis of permethylated glycosphingolipids is performed in the linear or reflector positive ion mode using α -dihydroxybenzoic acid (DHBA, 20mg/mL solution in 50%methanol: water) as a matrix. The procedures are described in detail below.

1. Prepare Matrix solution; 20mg α -dihydroxybenzoic acid (DHBA) in 1mL of 50% methanol: water in a 1.5 mL Microcentrifuge tube

2. Dissolve the permethylated sample with 20~ μL (Adjust the volume depends on your sample amount) of methanol using a clean glass micro syringe (rinsed with Methanol).
3. Take 1~2 μL of sample solution using a clean glass micro syringe (rinsed with Methanol), and then transfer into the bottom of a small vial.
4. Wait until the solution volume becomes approximately 1 μL (~3 min).
5. Take 1 μL of matrix solution using a micro pipet and then add into the sample solution.
6. Immediately mix them by pipetting the solution several times.
7. Take 1 μL of mixture and then put on a well of a MALDI plate. Immediately, take another 1ul of mixture then put on another well.
8. Wait until the mixtures on the plate dry.





Note:

- The matrix solution is preservable at -20°C . Defrost the matrix solution before use.
- A good crystallized sample is necessary to get a strong signal. However, making a good crystal is not always possible, hence, we recommend making at least two crystals per sample.

NSI-MSn analysis

Mass analysis by NSI-MS are performed by direct infusion of permethylated glycosphingolipids dissolved in infusion buffer the instrument at a constant flow rate of 0.4~1 $\mu\text{L}/\text{min}$ via a syringe pump. The procedures are descibed in detail below.

1. Prepare infusion buffer solution: Methanol/2-Propanol/1-Propanol/1 mM NaOH (16:3:3:2 by volume)
2. Dissolve the permethylated GSLs with 20~ μL (Adjust the volume depends on your sample amount) of methanol using a clean glass micro syringe (rinsed with methanol).
3. Take 10 μL of sample solution using a clean glass micro syringe (rinsed with methanol), and then transfer into a small vial (rinsed with methanol).
4. Add known concentration of heavy MeI permethylated Dp4 as a reference standard if quantification is desired (Optional).
5. Take 40~65 μL of the infusion buffer and then add to the sample solution.
6. Screw the cap, mix them by vortexing, and briefly centrifuge the vial.



Copyright:



Kazuhiro Aoki, 2018

© 2018 by Kazuhiro Aoki. This work is made available under the Creative Commons Attribution-NonCommercial-ShareAlike 4.0 International License: <http://creativecommons.org/licenses/by-nc-sa/4.0/>